

YJ-418 Veterinary Waterproof Scanner



- Image smoothing, filtering, edge enhancement, frame correlation, linear interpolation and other processing, the image is crisp and clear, high resolution.
- Measured for distance, perimeter, area, obstetric and etc.
- Individually adjustable near field, far field gain (5 stage adjustable).
- Single-B display, continuous real-time status in the field of vision to move the function.
- With display playback, gamma correction, histogram and other image processing functions.
- With built-in memory, and has an external USB interface, the built-in memory data transfer to PC.
- Light weight and convenient.

Standard configuration:

- Main Unit
- PS157 AC adapter
- CB2280 battery

- 6.5MHz rectal probe
- powder wire
- 250ml couplant
- Manual/ technical instructions
- bag with strap

Options:

- 7.5MHz linear probe
- 5.0MHz micro-convex probe
- 3.5MHz convex probe
- Video printer
- Detect Depth(mm):≥80
- Blind zone(mm):≤3
- Monitor size:7 inch
- Display mode B, B+B, B+M, M, 4B
- Image gray scale 256
- Scan range convex array
- Scan lines 512 lines/ frame
- Frame rate 30 frame/second
- Digital scan converter 512x512x8bits
- Body marks ≥15
- Multi-frequency 2MHz-10MHz
- Memory capacity ≥2GB
- Video output PAL
- Size 216mm(L) x 147mm(W) x 27mm(H)
- Net weight 900G
- Battery 3000Mah/7.4V
- Battery continuous work ≥3hours

Primary Competitive Advantages:

- Brand-name Parts
- Packaging
- Quality Approvals
- Origin
- Green Product
- Price
- Reputation
- Distributorships Offered

FORERMED®
www.yongjiuimed.com

- Guarantee/Warranty
- Product Features
- Service
- International Approvals
- Product Performance
- Small Orders Accepted
- Experienced Staff
- Prompt Delivery

Detail picture and Packing box:





